

KERN 347-117

Individual weight, OIML Class M1, 1 kg, Slot weight (Type B), stainless steel fine turned



Article codes, product groups & Status information

Model article-code	347-117
Model-Serie	347-117
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761366770
Planned start date	2024-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 kg
Tolerance (OIML)	50 mg

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	1 kg
OIML class	M1
OIML conform	Yes
Tolerance (OIML)	50 mg
Density (kg/m ³)	8.000
Volume (cm ³)	119,0476

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	118×118×17,1 mm
Gross weight	1,00 kg
Net weight	1,00 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-641
-------------------	---------

Functions

KERN 347-117

Individual weight, OIML Class M1, 1 kg, Slot weight (Type B), stainless steel fine turned



Standard



Option

