

PRODUCT BROCHURE

KERN 307-15

Date: 09.03.2026

ENGLISH



Individual weight, OIML Class E1, 20 kg, Knob, stainless steel polished

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Individual weight, OIML Class E1, 20 kg, Knob, Stainless steel polished
- Note: For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	307-15
Model-Serie	307
Assortment	KERN
Product Type	Test weights
Product group	OIML E1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761002791
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	20 kg
Tolerance (OIML)	10 mg
Case aluminium	317-150-600
Case wood	317-150-100

Technical data - Construction

Material	stainless steel polished
----------	--------------------------

Technical data - Measuring system

Nominal value	20 kg
OIML class	E1
OIML conform	yes
Tolerance (OIML)	10 mg
Density (kg/m ³)	8000
Volume (cm ³)	2500

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	129×129×227 mm
Gross weight	20 kg
Net weight	20 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	16 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	963-245
-------------------	---------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>317-150-100</u>	Wooden Box for Set of Weights KERN 317-150-100
<u>317-150-600</u>	Aluminium Protective Box KERN 317-150-600

SERVICES

Model	Description
<u>963-245</u>	DAkks Calibration KERN 963-245
<u>962-245R</u>	DAkks Recalibration KERN 962-245R
<u>921-245</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-245