

## PRODUCT BROCHURE

---

# OPTICS ODC-A8104

Date: 30.04.2025

ENGLISH



Eyepiece adapters for microscope cameras (0,37×/23,2 mm)

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## DESCRIPTION

Eyepiece adapter for microscope cameras ODC-A8104

## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                    |
|-----------------------|--------------------|
| Model article-code    | ODC-A8104          |
| Model-Serie           | ODC-AK             |
| Assortment            | Optics             |
| Product Type          | Microscopes        |
| Product group         | Microscope cameras |
| Product usage type    | Accessories/DL     |
| EAN code (Model)      | 4045761190825      |
| Customs tariff number | 90119000           |

### FUNCTIONS

Standard



ENGLISH

### Technical data - Packing & Shipping

|                                            |                |
|--------------------------------------------|----------------|
| Dimensions packaging (W×D×H)               | 90×90×90 mm    |
| Gross weight                               | 0.50 kg        |
| Net weight                                 | 0.20 kg        |
| Shipping method                            | Parcel service |
| Packaging component - by weight - plastics | 300 g          |
| Delivery time                              | 1 d            |

### Technical data - Environmental conditions

|                  |       |
|------------------|-------|
| Storage Temp Min | -5 °C |
| Stotage Temp Max | 40 °C |