

PRODUCT BROCHURE

OPTICS OBB-A1472

Date: 28.03.2025

ENGLISH



Binocular lens barrel

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

PRODUCT SHOWCASE



DESCRIPTION

Tube Binocular OBB-A1472

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OBB-A1472
Model-Serie	OBB-AT
Assortment	Optics
Product Type	Microscopes
Product group	Transmitted light microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761197473
Customs tariff number	90119000

Technical data - Primary

Tube	Binocular
------	-----------

Technical data - Construction

Tube type	Siedentopf
Tube angle	45°
Tube rotation	yes
Eye distance (max)	75 mm
Eye distance (min)	55 mm
Adjustable dioptre	both-sided

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	265×185×165 mm
------------------------------	----------------

Gross weight 1.6 kg

Net weight 1.2 kg

Shipping method Parcel service

Packaging component -
by weight - plastics 300 g

Delivery time 1 d

FUNCTIONS

Standard

