

## PRODUCT BROCHURE

# KERN SCD-4.0-DLS05

Date: 18.03.2025

ENGLISH



Flexible recording or transfer of measurements, particularly to Microsoft® Excel or Access

## CONTACT

**KERN & SOHN GmbH**

Ziegelei 1

72336 Balingen


Germany

Phone : +49 7433 9933-0


Fax : +49 7433 9933-149

Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)Web : [www.kern-sohn.com](http://www.kern-sohn.com)

## PRODUCT SHOWCASE



ENGLISH



## DESCRIPTION

- For operating systems Windows 10 (20H2) and newer, 11
- .NET Framework 4.8.1 required
- Full administration rights are required for the installation
- Supports balances and measuring devices with USB, RS-232, RS 485, Bluetooth, LAN or WiFi & 160;network (TCP/UDP/IP)
- Highly flexible formatting and design of the output (any order and positioning, formatting and rounding), particularly recording of date and time for every value transferred, if required
- Any number of devices/interfaces can be connected, as well as simultaneous and synchronised recording of several balances
- Key-activated or time-controlled interrogation of measurements or trigger of device functions, also for continuous recording
- The interface protocols for KERN balances are already predefined (standard configuration). Support of other balances is possible through flexible configuration
- Compatible interface cable included when you order a KERN balance at the same time

- Many different transfer and recording options:
  - Microsoft Excel/Access/Word in the foreground or background
  - Other Windows applications (through key simulation), e.g. shipping software or ERP system (SAP, Sage, etc.)
  - File recording (e.g. as text or CSV file)
  - Print out on text or label printer
- Screen output such as large display, line chart (drying curve for moisture analyzer), histogram, etc.
- : Benefits
  - Typing errors during manual transfer, for example, are avoided
  - Automatic recording, e.g. for long-term tests
  - GLP-compliant lab software
  - Can be used across different branches and applications
- Scope of supplies: Download link for 5 licenses

## SPECIFICATIONS

### Article codes, product groups & Status information

Model article-code	SCD-4.0-DLS05
Model-Serie	SCD
Assortment	KERN
Product Type	Software
Product group	BalanceConnection
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1488700
Predecessor Model	SCD-4.0S05
Customs tariff number	852349 10

### Technical data - Packing & Shipping

Delivery time	1 d
---------------	-----

## SERVICES

Model	Description
970-105	Warranty extension KERN (+2 years) 970-105