

KERN 346-07

Individual weight, OIML Class M1, 10 kg, Block, stainless steel



Article codes, product groups & Status information

Model article-code	346-07
Model-Serie	346
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761117136
Planned start date	2010-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	10 kg
Tolerance (OIML)	500 mg
Case aluminium	346-070-600

Technical data - Construction

Dimensions	193×96×117 mm
Material	Stainless steel

Technical data - Measuring system

Nominal value	10 kg
OIML class	M1
OIML conform	Yes
Tolerance (OIML)	500 mg
Density (kg/m ³)	7200
Volume (cm ³)	1389

Functions

Standard

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	193×96×117 mm
Gross weight	10 kg
Net weight	10 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-644
-------------------	---------

KERN 346-07

Individual weight, OIML Class M1, 10 kg, Block, stainless steel

