

KERN 326-01

Individual weight, OIML Class F1, 1 g, Eco Shape, stainless steel polished



Article codes, product groups & Status information

Model article-code	326-01
Model-Serie	326
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761223745
Predecessor Model	326-01-2018a; 327-61; 329-01
Successor Model	326-01-2025e
Planned start date	2019-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	1 g
Tolerance (OIML)	0,1 mg
Case aluminium	317-010-600
Case plastic	347-030-400
Case wood	317-010-100

Technical data - Construction

Material	Stainless steel polished
----------	--------------------------

Technical data - Measuring system

Nominal value	1 g
OIML class	F1
OIML conform	Yes

Functions

Tolerance (OIML)	0,1 mg
------------------	--------

Density (kg/m ³)	7950
------------------------------	------

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	6×6×6 mm
------------------------------	----------

Gross weight	0.05 kg
--------------	---------

Net weight	0.05 kg
------------	---------

Shipping method	Parcel service
-----------------	----------------

Packaging component - by weight - plastics	0 g
--	-----

Delivery time	1 d
---------------	-----

Optional Services

DAkKS Calibration	962-431
-------------------	---------

KERN 326-01

Individual weight, OIML Class F1, 1 g, Eco Shape, stainless steel polished



Standard

