


PRODUCT BROCHURE

OPTICS OBB-A1259

Date: 30.04.2025

ENGLISH



Infinity planachromatic objective lens (without cover glass) for a large working distance, 40 x / 0,65, (3,9 mm) (sprung)

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

Objective lens 40 x (Infinity Plan) OBB-A1259

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OBB-A1259
Model-Serie	OBB-AO
Assortment	Optics
Product Type	Microscopes
Product group	Metallurgical microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761167940
Customs tariff number	90119000

FUNCTIONS

Standard



ENGLISH

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	65×40×40 mm
Gross weight	0.60 kg
Net weight	0.30 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	300 g
Delivery time	1 d