

PRODUCT BROCHURE

KERN 367-07

Date: 10.03.2026

ENGLISH



Individual weight, OIML Class M3, 100 g, Knob, stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	367-07
Model-Serie	367-07
Assortment	KERN
Product Type	Test weights
Product group	OIML M3 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761257269
Predecessor Model	367-47
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	100 g
Tolerance (OIML)	50 mg
Case aluminium	317-070-600
Case plastic	347-070-400
Case wood	337-070-200

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	100 g
---------------	-------

OIML class M3

OIML conform yes

Tolerance (OIML) 50 mg

Density (kg/m³) 7950

Volume (cm³) 11,9048

Technical data - Packing & Shipping

Dimensions packaging (W×D×H) 22×22×37 mm

Gross weight 0.10 kg

Net weight 0.10 kg

Shipping method Parcel service

Packaging component - by weight - plastics 0 g

Delivery time 1 d

Technical data - Verification

Verification possible no

Optional Services

DAkkS Calibration 962-637

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>337-070-200</u>	Wooden Box for Set of Weights KERN 337-070-200
<u>347-070-400</u>	Plastic Weight Box KERN 347-070-400
<u>317-070-600</u>	Aluminium Protective Box KERN 317-070-600
<u>315-242</u>	Tweezers KERN 315-242

SERVICES

Model	Description
<u>962-637</u>	DAkkS Calibration KERN 962-637
<u>962-637R</u>	DAkkS Recalibration KERN 962-637R