

PRODUCT BROCHURE

KERN 323-044

Date: 09.03.2026

ENGLISH



Weight set, F1, 1 mg - 200 g, Knob, stainless steel polished, in plastic case

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



DESCRIPTION

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	323-044
Model-Serie	323
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main article
EAN code (Model)	404576 1080737
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 mg - 200 g
---------------	--------------

Case plastic	313-052-400
Weight set case	in plastic case

Technical data - Construction

Material	stainless steel polished
Is weight set	yes
Number of weights in set	23
Type of case	in plastic case

Technical data - Measuring system

Nominal value	1 mg - 200 g
OIML class	F1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	250×210×75 mm
Gross weight	0.95 kg
Net weight	0.95 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-403
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>962-403</u>	DAkkS Calibration KERN 962-403
<u>962-403R</u>	DAkkS Recalibration KERN 962-403R
<u>921-403</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-403

Model	Description
<u>920-704</u>	Marking (Last digits of serial number) KERN 920-704

SPARE PARTS

Model	Description
<u>313-052-400</u>	Plastic Weight Box KERN 313-052-400