

KERN 347-136

Individual weight, OIML Class M1, 5 kg, Hook weight, stainless steel fine turned



Article codes, product groups & Status information

Model article-code	347-136
Model-Serie	347-136
Assortment	Scales
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761255418
Predecessor Model	347-536
Planned start date	2022-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	5 kg
Tolerance (OIML)	250 mg
Case plastic	347-140-400

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	5 kg
OIML class	M1
OIML conform	Yes
Tolerance (OIML)	250 mg
Density (kg/m ³)	8.000

Functions

Standard

Volume (cm³) 595,2381

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	80×80×130 mm
Gross weight	5,0 kg
Net weight	5,0 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkKS Calibration	962-643
-------------------	---------

KERN 347-136

Individual weight, OIML Class M1, 5 kg, Hook weight, stainless steel fine turned



Option

