

PRODUCT BROCHURE

KERN 367-11

Date: 10.03.2026

ENGLISH



Individual weight, OIML Class M3, 1 kg, Knob, stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	367-11
Model-Serie	367-11
Assortment	KERN
Product Type	Test weights
Product group	OIML M3 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761257290
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 kg
Tolerance (OIML)	500 mg
Case aluminium	317-110-600
Case plastic	347-110-400
Case wood	337-110-200

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	1 kg
OIML class	M3

OIML conform	yes
Tolerance (OIML)	500 mg
Density (kg/m ³)	7950
Volume (cm ³)	119,0476

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	48×48×80 mm
Gross weight	1.00 kg
Net weight	1.00 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	no
-----------------------	----

Optional Services

DAkkS Calibration	962-641
-------------------	---------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>317-280</u>	Gloves KERN 317-280
<u>337-110-200</u>	Wooden Box for Set of Weights KERN 337-110-200
<u>347-110-400</u>	Plastic Weight Box KERN 347-110-400
<u>317-110-600</u>	Aluminium Protective Box KERN 317-110-600

SERVICES

Model	Description
<u>962-641</u>	DAkkS Calibration KERN 962-641
<u>962-641R</u>	DAkkS Recalibration KERN 962-641R