

PRODUCT BROCHURE

KERN 347-115

Date: 19.09.2025

ENGLISH



Individual weight, OIML Class M1, 1 kg, Slot weight (Type A), stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Individual weight, OIML Class M1, 1 kg, Slot weight (Type A), stainless steelfine turned (OIML)

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	347-115
Model-Serie	347-115
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 125636 1
Predecessor Model	347-515
Customs tariff number	84239090

Technical data - Primary

Nominal value	1 kg
Tolerance (OIML)	50 mg
Case plastic	347-130-400

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	1 kg
OIML class	M1

OIML conform	yes
Tolerance (OIML)	50 mg
Density (kg/m ³)	8000
Volume (cm ³)	119,0476

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	80×80×28 mm
Gross weight	1.00 kg
Net weight	1.00 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-641
-------------------	---------

FUNCTIONS

Standard



ACCESSORIES

Model	Description
317-280	Gloves KERN 317-280
347-130-400	Plastic Weight Box KERN 347-130-400
347-075-100	Supporting Bar KERN 347-075-100
347-095-100	Supporting Bar KERN 347-095-100
347-115-100	Supporting Bar KERN 347-115-100

SERVICES

Model	Description
962-641	DAkkS Calibration KERN 962-641
962-641R	DAkkS Recalibration KERN 962-641R
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351

RELATED ARTICLES DATA

Model name	Nominal value	Tolerance (OIML)	Case plastic
347-115	1 kg	50 mg	347-130-400