

PRODUCT BROCHURE

KERN 303-064

Date: 10.03.2026

ENGLISH



Weight set, E1, 1 mg - 1 kg, Knob, stainless steel polished, in plastic case

CONTACT


KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



DESCRIPTION

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	303-064
Model-Serie	303
Assortment	KERN
Product Type	Test weights
Product group	OIML E1 Test weights
Product usage type	Main article
EAN code (Model)	4045761291614
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 mg - 1 kg
---------------	-------------

Case plastic	313-082-400
Weight set case	in plastic case

Technical data - Construction

Material	stainless steel polished
Is weight set	yes
Number of weights in set	25
Type of case	in plastic case

Technical data - Measuring system

Nominal value	1 mg - 1 kg
OIML class	E1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	350×320×180 mm
Gross weight	3.8 kg
Net weight	3.8 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	963-205
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>963-205</u>	DAkkS Calibration KERN 963-205
<u>962-205R</u>	DAkkS Recalibration KERN 962-205R
<u>920-714</u>	Marking (Last digits of serial number) KERN 920-714

SPARE PARTS

Model	Description
<u>313-082-400</u>	Plastic Weight Box KERN 313-082-400
<u>315-242</u>	Tweezers KERN 315-242
<u>318-270</u>	Dusting Brush KERN 318-270