

PRODUCT BROCHURE

KERN 323-074

Date: 09.03.2026

ENGLISH



Weight set, F1, 1 mg - 2 kg, Knob, stainless steel polished, in plastic case


CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information



ENGLISH

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

Model article-code	323-074
Model-Serie	323
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main article
EAN code (Model)	404576 1080706
Customs tariff number	84239090
Country of Origin	DE

Nominal value	1 mg - 2 kg
---------------	-------------

Case plastic	313-082-400
Weight set case	in plastic case

Material	stainless steel polished
Is weight set	yes
Number of weights in set	27
Type of case	in plastic case

Nominal value	1 mg - 2 kg
OIML class	F1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	350×320×180 mm
Gross weight	8 kg
Net weight	8 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-406
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>962-406</u>	DAkkS Calibration KERN 962-406
<u>962-406R</u>	DAkkS Recalibration KERN 962-406R
<u>921-406</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-406

Model	Description
<u>920-707</u>	Marking (Last digits of serial number) KERN 920-707

SPARE PARTS

Model	Description
<u>313-082-400</u>	Plastic Weight Box KERN 313-082-400