

# KERN 326-03

Individual weight, OIML Class F1, 5 g, Eco Shape, stainless steel polished



## Article codes, product groups & Status information

|                       |                                |
|-----------------------|--------------------------------|
| Model article-code    | 326-03                         |
| Model-Serie           | 326                            |
| Assortment            | KERN                           |
| Product Type          | Test weights                   |
| Product group         | OIML F1 Test weights           |
| Product usage type    | Main item;Accessories/Services |
| EAN code (Model)      | 4045761223721                  |
| Predecessor Model     | 326-03-2018a; 327-63; 329-03   |
| Successor Model       | 326-03-2025e                   |
| Planned start date    | 2019-1HJ                       |
| Customs tariff number | 84239090                       |

## Technical data - Primary

|                  |             |
|------------------|-------------|
| Nominal value    | 5 g         |
| Tolerance (OIML) | 0,16 mg     |
| Case aluminium   | 317-030-600 |
| Case plastic     | 347-030-400 |
| Case wood        | 317-030-100 |

## Technical data - Construction

|          |                          |
|----------|--------------------------|
| Material | Stainless steel polished |
|----------|--------------------------|

## Technical data - Measuring system

|               |     |
|---------------|-----|
| Nominal value | 5 g |
| OIML class    | F1  |
| OIML conform  | Yes |

## Functions

Tolerance (OIML) 0,16 mg

Density (kg/m<sup>3</sup>) 7950

## Technical data - Packing & Shipping

Dimensions packaging (W×D×H) 8×8×13 mm

Gross weight 0.05 kg

Net weight 0.05 kg

Shipping method Parcel service

Packaging component - by weight - plastics 0 g

Delivery time 1 d

## Optional Services

DAkKS Calibration 962-433

# KERN 326-03

Individual weight, OIML Class F1, 5 g, Eco Shape, stainless steel polished



Standard

