

PRODUCT BROCHURE

KERN 303-024

Date: 10.03.2026

ENGLISH



Weight set, E1, 1 mg - 50 g, Knob, stainless steel polished, in plastic case


CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information



ENGLISH

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

Article codes, product groups & Status information

Case plastic	313-052-400
Weight set case	in plastic case

Material	stainless steel polished
Is weight set	yes
Number of weights in set	20
Type of case	in plastic case

Nominal value	1 mg - 50 g
OIML class	E1
OIML conform	yes

Nominal value	1 mg - 50 g
---------------	-------------

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	250×210×75 mm
Gross weight	0.80 kg
Net weight	0.80 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	963-201
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>963-201</u>	DAkkS Calibration KERN 963-201
<u>962-201R</u>	DAkkS Recalibration KERN 962-201R
<u>920-710</u>	Marking (Last digits of serial number) KERN 920-710

SPARE PARTS

Model	Description
<u>313-052-400</u>	Plastic Weight Box KERN 313-052-400
<u>315-242</u>	Tweezers KERN 315-242
<u>318-270</u>	Dusting Brush KERN 318-270