

## PRODUCT BROCHURE

---

# KERN 317-02-2020a

Date: 15.03.2025

ENGLISH



Individual weight, OIML Class E2, 2 g, Knob, stainless steel polished

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)

## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                          |
|-----------------------|--------------------------|
| Model article-code    | 317-02-2020a             |
| Model-Serie           | 317                      |
| Assortment            | KERN                     |
| Product Type          | Test weights             |
| Product group         | OIML E2 Test weights     |
| Product usage type    | Trade parts list element |
| EAN code (Model)      | 404576 1003545           |
| Successor Model       | 317-02; 317-02-2022i     |
| Customs tariff number | 84239090                 |

### Technical data - Primary

|                  |             |
|------------------|-------------|
| Nominal value    | 2 g         |
| Tolerance (OIML) | 0,04 mg     |
| Case aluminium   | 317-020-600 |
| Case plastic     | 317-020-400 |
| Case wood        | 317-020-100 |

### Technical data - Construction

|          |                          |
|----------|--------------------------|
| Material | stainless steel polished |
|----------|--------------------------|

### Technical data - Measuring system

|                              |      |
|------------------------------|------|
| OIML class                   | E2   |
| OIML conform                 | yes  |
| Density (kg/m <sup>3</sup> ) | 7950 |
| Volume (cm <sup>3</sup> )    | 0,25 |

### Technical data - Packing & Shipping

|                                            |                |
|--------------------------------------------|----------------|
| Dimensions packaging (W×D×H)               | 6×6×10 mm      |
| Gross weight                               | 0.05 kg        |
| Net weight                                 | 0.05 kg        |
| Shipping method                            | Parcel service |
| Packaging component - by weight - plastics | 0 g            |
| Delivery time                              | 1 d            |

### Technical data - Verification

|                       |     |
|-----------------------|-----|
| Verification possible | yes |
|-----------------------|-----|

### Optional Services

|                   |         |
|-------------------|---------|
| DAkkS Calibration | 962-332 |
|-------------------|---------|

## FUNCTIONS

|          |
|----------|
| Standard |
|----------|



## ACCESSORIES

| Model       | Description                                    |
|-------------|------------------------------------------------|
| 317-020-100 | Wooden Box for Set of Weights KERN 317-020-100 |
| 317-020-400 | Plastic Weight Box KERN 317-020-400            |
| 317-020-600 | Aluminium Protective Box KERN 317-020-600      |
| 315-242     | Tweezers KERN 315-242                          |

## SERVICES

| Model    | Description                                                        |
|----------|--------------------------------------------------------------------|
| 952-332  | Verification KERN 952-332                                          |
| 962-332  | DAkkS Calibration KERN 962-332                                     |
| 962-332R | DAkkS Recalibration KERN 962-332R                                  |
| 920-33   | Marking by etching KERN 920-33                                     |
| 921-332  | Guaranteed Plus-deviation ( $\pm 1/3$ OIML) KERN 921-332           |
| 920-331  | Marking by etching (last digits of the serial number) KERN 920-331 |