

KERN CFS 6K0.1

Professional, high-resolution counting scale with 100 item memories and second balance interface, counting resolution up to 75,000 points



Article codes, product groups & Status information

Model article-code	CFS 6K0.1
Model-Serie	CFS
Assortment	KERN
Product Type	Industrial scales
Product group	Counting scales
Product usage type	Main article
EAN code (Model)	4045761124332
Planned start date	2011-1HJ
Customs tariff number	84238129
CE Label	Yes

Technical data - Primary

Weighing capacity	6 kg
Readability	0.0001 kg
Counting resolution	60000
Minimum piece weight	1 g
Net weight	3.2 kg
Platform size	295×225 mm

Technical data - Construction

Casing material	Plastic
Dimensions of housing	315×350×100 mm
Dimensions device	320×350×125 mm
Material weighing plate	Stainless steel
Bubble Level	Yes

Revolving screw feet	Yes
Weighing surface (d)	295×225 mm
Dimensions of weighing plate (W×D)	295×225 mm
IP protection load cell	IP65

Technical data - Display

Display height	20 mm
Backlight	Yes

Technical data - Power supply

Mains adapter/adaptor included	CH;EURO;UK;US
Charging time	14 h
Input Voltage	220 V - 240 V AC, 50 Hz
Operating time (Backlight off)	70 h
Rechargeable battery	Rchrg. battery optional

Technical data - Measuring system

Calibration method	Adjusting with external weight
Linearity	± 0.0002 kg
Readability	0.0001 kg
Recommended adjusting weight	1 kg (F2); 5 kg (F2)
Repeatability	0.0001 kg
Resolution	60.000
Stabilisation time under laboratory conditions	2 s
Tare range	6 kg
Warm up time	120 min

KERN CFS 6K0.1

Professional, high-resolution counting scale with 100 item memories and second balance interface, counting resolution up to 75,000 points



Weighing system	Strain gauge
-----------------	--------------

Weighing units	Kg;lb
----------------	-------

Technical data - Environmental conditions

Maximum humidity	80 %
------------------	------

Maximum operating temperature	40 °C
-------------------------------	-------

Minimum ambient temperature	0 °C
-----------------------------	------

Minimum humidity	0 %
------------------	-----

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	455×380×215 mm
------------------------------	----------------

Gross weight	4.4 kg
--------------	--------

Net weight	3.2 kg
------------	--------

Shipping method	Parcel service
-----------------	----------------

Packaging component - by weight - cardboard	537 g
---	-------

Packaging component - by weight - styrofoam	179 g
---	-------

Delivery time	1 d
---------------	-----

Technical data - Counting data

Counting resolution	60000
---------------------	-------

Minimum piece weight at piece counting (Laboratory)	100 mg
---	--------

Minimum piece weight	1 g
----------------------	-----

Optional Services

DAkkS Calibration	963-128
-------------------	---------

Functions

Standard



Option

