

PRODUCT BROCHURE

KERN 347-035

Date: 28.03.2025

ENGLISH



Individual weight, OIML Class M1, 5 g, Slot weight (Type A), stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

- Individual weight, OIML Class M1, 5 g, Slot weight (Type A), stainless steelfine turned (OIML)

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	347-035
Model-Serie	347-035
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1256279
Predecessor Model	347-435
Customs tariff number	84239090

Technical data - Primary

Nominal value	5 g
Tolerance (OIML)	1,6 mg
Case plastic	347-030-400

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	5 g
OIML class	M1

OIML conform	yes
Tolerance (OIML)	1,6 mg
Density (kg/m ³)	8000
Volume (cm ³)	0,5952

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	18×18×4 mm
Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-633
-------------------	---------

FUNCTIONS

Standard



ACCESSORIES

Model	Description
347-030-400	Plastic Weight Box KERN 347-030-400
347-075-100	Supporting Bar KERN 347-075-100
347-095-100	Supporting Bar KERN 347-095-100
347-115-100	Supporting Bar KERN 347-115-100
315-242	Tweezers KERN 315-242

SERVICES

Model	Description
962-633	DAkkS Calibration KERN 962-633
962-633R	DAkkS Recalibration KERN 962-633R
920-30	Marking with stroke number KERN 920-30
920-32	Marking by etching KERN 920-32
920-301	Marking with stroke number (Last digits of serial number) KERN 920-301
920-321	Marking by etching (last digits of the serial number) KERN 920-321

RELATED ARTICLES DATA

Model name	Nominal value	Tolerance (OIML)	Case plastic
347-035	5 g	1,6 mg	347-030-400