


PRODUCT BROCHURE

SAUTER AD 0032

Date: 30.04.2025

ENGLISH



Screw Tension Clamp without jaws

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	AD 0032
Model-Serie	AD
Assortment	Sauter
Product Type	Force measurement
Product group	Clamps & clamping tools
Product usage type	Trading parts list element;Accessories/Se rvices
EAN code (Model)	4045761174467
Customs tariff number	90278990

Dimensions packaging (W×D×H)	160×75×30 mm
---------------------------------	--------------

Gross weight	0.70 kg
--------------	---------

Net weight	0.70 kg
------------	---------

Shipping method	Parcel service
-----------------	----------------

UN number	No dangerous goods
-----------	--------------------

Packaging component - by weight - plastics	0 g
---	-----

Delivery time	14 d
---------------	------

ENGLISH

Technical data - Primary

Max. Force	1 kN
Max. Clamp width	100 mm

Technical data - Construction

Material	aluminium
Dimension housing (W×D×H)	160×75×30 mm

Technical data - Environmental conditions

Maximum storage temperature	40 °C
Minimum storage temperature	-10 °C

Technical data - Packing & Shipping

ACCESSORIES

Model	Description
AD 0034	Jaws KERN AD 0034
AD 0035	Jaws KERN AD 0035
AD 0036	Jaws KERN AD 0036
AD 0037	Jaws KERN AD 0037
AD 0038	Jaws KERN AD 0038
AD 0039	Jaws KERN AD 0039
AD 0040	Jaws KERN AD 0040
AD 0041	Jaws KERN AD 0041
AD 0042	Jaws KERN AD 0042
AD 0043	Jaws KERN AD 0043
AD 0044	Jaws KERN AD 0044
AD 0045	Jaws KERN AD 0045
AD 0046	Jaws KERN AD 0046