


KERN 318-07

Milligram weights, E2, 100 mg, Platelet, stainless steel polished



Article codes, product groups & Status information

Model article-code	318-07
Model-Serie	318
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761240469
Predecessor Model	318-7-2.019a
Planned start date	2020-1HJ
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	100 mg
Tolerance (OIML)	0,016 mg
Case aluminium	317-009-600
Case plastic	347-009-400
Case wood	338-090-200

Technical data - Construction

Dimensions (W×D×H)	8,5×8,5×4 mm
Material	Stainless steel polished

Technical data - Measuring system

Nominal value	100 mg
OIML class	E2

Functions

Standard

OIML conform	Yes
Tolerance (OIML)	0,016 mg
Density (kg/m ³)	7.960
Volume (cm ³)	0,0116

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	8,5×8,5×4 mm
Gross weight	0,05 kg
Net weight	0,05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-357
-------------------	---------

KERN 318-07

Milligram weights, E2, 100 mg, Platelet, stainless steel polished



Option

