

PRODUCT BROCHURE

KERN 315-103-188

Date: 10.03.2026

ENGLISH



Safety Set 5 g (E2); 200 mg (E2)

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



DESCRIPTION

- These “KERN Safety Sets” each contain one test weight for testing the sensitivity, i.e. the correct adjustment of your balance as well as a small test weight to test the lower end of the weighing range, the so-called minimum sample weight
- As an option, the “KERN Safety Set” has space for another test weight, for testing your balance at a weight which is relevant for you
- Useful accessories which have been selected to suit that particular “KERN Safety Set”, such as, for example, special gloves, tweezers, weight grips, brushes, etc. will assist you in handling your test weights properly
- Stored in the practical protective case next to your balance, you can check and ensure the high precision of your balance at any time

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	315-103-188
Model-Serie	315
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights

Product usage type	Main item;Accessories/Services
--------------------	--------------------------------

EAN code (Model)	4045761223356
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	200 mg & 5 g
---------------	--------------

Case aluminium	zet315-000-000
Weight set case	in aluminium case

Technical data - Construction

Material	stainless steel
Is weight set	yes
Number of weights in set	2
Type of case	in aluminium case

Technical data - Measuring system

Nominal value	200 mg & 5 g
OIML class	E2
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	150×115×70 mm
Gross weight	1.00 kg
Net weight	0.50 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	500 g
Delivery time	1 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-300-300
-------------------	-------------

FUNCTIONS

Standard



SERVICES

Model	Description
<u>921-334</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-334
<u>962-300-300</u>	DAkkS Calibration KERN 962-300-300

SPARE PARTS

Model	Description
<u>315-242</u>	Tweezers KERN 315-242