

## PRODUCT BROCHURE

---

# KERN 343-06

Date: 09.03.2026

ENGLISH



Weight set, M1, 1 mg - 1 kg, Knob, stainless steel fine turned, in wooden case

---

### CONTACT


KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## PRODUCT SHOWCASE



## DESCRIPTION

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                       |
|-----------------------|-----------------------|
| Model article-code    | 343-06                |
| Model-Serie           | 343                   |
| Assortment            | KERN                  |
| Product Type          | Test weights          |
| Product group         | OIML M 1 Test weights |
| Product usage type    | Main article          |
| EAN code (Model)      | 4045761007543         |
| Customs tariff number | 84239090              |
| Country of Origin     | DE                    |

## Technical data - Primary

|               |             |
|---------------|-------------|
| Nominal value | 1 mg - 1 kg |
|---------------|-------------|

|                 |                |
|-----------------|----------------|
| Case wood       | 353-460-200    |
| Weight set case | in wooden case |

## Technical data - Construction

|                          |                             |
|--------------------------|-----------------------------|
| Material                 | stainless steel fine turned |
| Is weight set            | yes                         |
| Number of weights in set | 25                          |
| Type of case             | in wooden case              |

## Technical data - Measuring system

|               |             |
|---------------|-------------|
| Nominal value | 1 mg - 1 kg |
| OIML class    | M1          |
| OIML conform  | yes         |

## Technical data - Packing & Shipping

|                                            |                |
|--------------------------------------------|----------------|
| Dimensions packaging (W×D×H)               | 230×180×112 mm |
| Gross weight                               | 4.4 kg         |
| Net weight                                 | 4.4 kg         |
| Shipping method                            | Parcel service |
| Packaging component - by weight - plastics | 0 g            |
| Delivery time                              | 1 d            |

## Optional Services

|                   |         |
|-------------------|---------|
| DAkkS Calibration | 962-605 |
|-------------------|---------|

## FUNCTIONS

|          |
|----------|
| Standard |
|----------|



|        |
|--------|
| Option |
|--------|



## SERVICES

| Model                           | Description                                         |
|---------------------------------|-----------------------------------------------------|
| <a href="#"><u>962-605</u></a>  | DAkkS Calibration KERN 962-605                      |
| <a href="#"><u>962-605R</u></a> | DAkkS Recalibration KERN 962-605R                   |
| <a href="#"><u>920-706</u></a>  | Marking (Last digits of serial number) KERN 920-706 |

## SPARE PARTS

| Model                              | Description                     |
|------------------------------------|---------------------------------|
| <a href="#"><u>353-460-200</u></a> | Box of Weights KERN 353-460-200 |
| <a href="#"><u>315-242</u></a>     | Tweezers KERN 315-242           |
| <a href="#"><u>318-270</u></a>     | Dusting Brush KERN 318-270      |