

# KERN 347-036-N00

Individual weight (Newton), 5 cN, Hook weight, stainless steel fine turned



## Article codes, product groups & Status information

|                       |                                |
|-----------------------|--------------------------------|
| Model article-code    | 347-036-N00                    |
| Model-Serie           | 347-036-N00                    |
| Assortment            | KERN                           |
| Product Type          | Test weights                   |
| Product group         | OIML M1 Test weights           |
| Product usage type    | Main item;Accessories/Services |
| EAN code (Model)      | 4045761262492                  |
| Planned start date    | 2022-1HJ                       |
| Customs tariff number | 84239090                       |
| Country of Origin     | DE                             |

## Technical data - Primary

|                  |             |
|------------------|-------------|
| Tolerance (OIML) | 1,6 mg      |
| Case plastic     | 347-030-400 |

## Technical data - Construction

|          |                             |
|----------|-----------------------------|
| Material | Stainless steel fine turned |
|----------|-----------------------------|

## Technical data - Measuring system

|                              |        |
|------------------------------|--------|
| OIML class                   | M1     |
| OIML conform                 | No     |
| Tolerance (OIML)             | 1,6 mg |
| Density (kg/m <sup>3</sup> ) | 7.950  |
| Volume (cm <sup>3</sup> )    | 0,5952 |

## Technical data - Packing & Shipping

|                              |           |
|------------------------------|-----------|
| Dimensions packaging (W×D×H) | 8×8×22 mm |
|------------------------------|-----------|

|                                            |                |
|--------------------------------------------|----------------|
| Gross weight                               | 0,05 kg        |
| Net weight                                 | 0,05 kg        |
| Shipping method                            | Parcel service |
| Packaging component - by weight - plastics | 0 g            |
| Delivery time                              | 6 d            |

## Optional Services

|                   |          |
|-------------------|----------|
| DAkkS Calibration | 962-633N |
|-------------------|----------|

## Functions

# KERN 347-036-N00

Individual weight (Newton), 5 cN, Hook weight, stainless steel fine turned



## Standard



## Option

