

KERN 572-33

All-rounder, e.g. as precision balance in the laboratory or in harsh industrial applications, ideal for the diverse possibilities of Industry 4.0 applications

KERN®

Article codes, product groups & Status information

Model article-code	572-33
Model-Serie	572
Assortment	KERN
Product Type	Laboratory balances
Product group	Precision balances
Product usage type	Main article
EAN code (Model)	4045761356856
Predecessor Model	572-33-2021a
Planned start date	2022-2HJ
Customs tariff number	90160010
CE Label	Yes

Technical data - Primary

Weighing capacity	1600 g
Readability	0.01 g
Repeatability	0.01 g
Linearity	± 0.03 g
Minimum piece weight	100 mg
Weighing surface (d)	150 mm

Technical data - Construction

Dimensions of housing	180×310×85 mm
Bubble Level	Yes
Material weighing plate	Stainless steel
Revolving screw feet	Yes

Casing material	Cast aluminium
Weighing surface (d)	150 mm
Dimensions of weighing plate (diameter)	150 mm

Technical data - Display

Display height	21 mm
Backlight	Yes

Technical data - Power supply

Mains adapter/adaptor included	CH;EURO;UK;US
Charging time	8 h
Input Voltage	100 V - 240 V AC, 50 / 60 Hz
Operating time (Backlight off)	48 h
Rechargeable battery	Rchrg. battery optional (factory)

Technical data - Measuring system

Calibration method	Adjusting with external weight
Linearity	± 0.03 g
Minimum weight (USP)	20 g
Readability	0.01 g
Recommended adjusting weight	1 kg (F1); 500 g (F1)
Repeatability	0.01 g
Stabilisation time under laboratory conditions	3 s
Warm up time	120 min

Technical data - Environmental conditions

KERN 572-33

All-rounder, e.g. as precision balance in the laboratory or in harsh industrial applications, ideal for the diverse possibilities of Industry 4.0 applications



Maximum humidity 80 %

Maximum operating temperature 40 °C

Minimum ambient temperature -10 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H) 400×340×220 mm

Gross weight 4.4 kg

Net weight 2.4 kg

Shipping method Parcel service

Packaging component - by weight - cardboard 1050 g

Packaging component - by weight - plastics 1050 g

Delivery time 1 d

Technical data - Counting data

Counting resolution 160000

Minimum piece weight at piece counting (Laboratory) 10 mg

Minimum piece weight 100 mg

Optional Services

DAkKS Calibration 963-127

Functions

Standard



Option



Factory

