

PRODUCT BROCHURE

KERN 347-055-N00

Date: 28.03.2025

ENGLISH



Individual weight (Newton), 20 cN, Slot weight (Type A), stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	347-055-N00
Model-Serie	347-055-N00
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1277 175
Customs tariff number	84239090

Technical data - Primary

Tolerance (OIML)	2,5 mg
Case plastic	347-080-400

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

OIML class	M1
OIML conform	no
Tolerance (OIML)	2,5 mg
Density (kg/m ³)	7950
Volume (cm ³)	2,381

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	27×27×6 mm
------------------------------	------------

Gross weight	0.05 kg
--------------	---------

Net weight	0.05 kg
------------	---------

Shipping method	Parcel service
-----------------	----------------

Packaging component - by weight - plastics	0 g
--	-----

Delivery time	1 d
---------------	-----

Optional Services

DAkkS Calibration	962-635N
-------------------	----------

FUNCTIONS

Standard



ACCESSORIES

Model	Description
347-080-400	Plastic Weight Box KERN 347-080-400
347-075-100-N00	Supporting Bar KERN 347-075-100-N00
347-095-100-N00	Supporting Bar KERN 347-095-100-N00
347-115-100-N00	Supporting Bar KERN 347-115-100-N00
315-242	Tweezers KERN 315-242

SERVICES

Model	Description
962-635RN	DAkkS Recalibration KERN 962-635RN
920-32	Marking by etching KERN 920-32
962-635N	DAkkS Calibration KERN 962-635N
920-301	Marking with stroke number (Last digits of serial number) KERN 920-301
920-321	Marking by etching (last digits of the serial number) KERN 920-321

RELATED ARTICLES DATA

Model name	Tolerance (OIML)	Case plastic
347-055-N00	2,5 mg	347-080-400