



Article codes, product groups & Status information

Model article-code	TVP
Model-Serie	TVP
Assortment	Testing & Measuring Devices
Product Type	Force measurement
Product group	Test stands (mechanic)
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761139756
Predecessor Model	TVP-2.019a
Planned start date	2012-1HJ
Customs tariff number	90248000
Country of Origin	CN

Technical data - Construction

Dimension housing (W×D×H)	150×233×420 mm
Mounting threads	5x threaded hole M6
Measuring room height	315 mm

Technical data - Environmental conditions

Maximum storage temperature	40 °C
Minimum storage temperature	-10 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	575×238×280 mm
Gross weight	12 kg
Net weight	11 kg

Shipping method

UN number	Parcel service
Packaging component - by weight - cardboard	No dangerous goods
Packaging component - by weight - styrofoam	750 g
Delivery time	750 g
	1 d

KERN TVP

Manual test stand for compressive force measurements



Functions

Standard

