

KERN OZB-A7101EU

Professional illumination units guarantee outstanding, uniform and strong illumination



Article codes, product groups & Status information

Model article-code	OZB-A7101EU
Model-Serie	OZB-IR
Assortment	Optics
Product Type	Microscopes
Product group	External Illumination
Product usage type	Main article
EAN code (Model)	4045761361584
Predecessor Model	OZB-A7101
Planned start date	2018-1HJ
Customs tariff number	85395100

Technical data - Primary

Illumination type	LED
Illumination intensity	4.5W

Technical data - Construction

Color tempertaure	6500-7000 K
Dimmable	Yes
Ring diameter	6.2 cm
Polarization filter	Yes
Dimensions	99×90×33 mm

Technical data - Power supply

Protection Class	Class III - Low voltage
------------------	-------------------------

Technical data - Environmental conditions

Storage Temp Min	0 °C
Stotage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	250×150×80 mm
Gross weight	1.00 kg
Net weight	0.95 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	62.5 g
Packaging component - by weight - plastics	0 g
Packaging component - by weight - styrofoam	0 g
Packaging component - by weight - foam	0 g
Delivery time	1 d

Technical data - Verification

WEEE	Yes
------	-----

Functions

Standard

KERN OZB-A7101EU

Professional illumination units guarantee outstanding, uniform and strong illumination

