

PRODUCT BROCHURE

OPTICS ODC 825

Date: 15.12.2025

ENGLISH



Specialists in microscopy for measurement, counting, documentation, archiving and image processing

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	ODC 825
Model-Serie	ODC-82
Assortment	Optics
Product Type	Microscopes
Product group	Microscope cameras
Product usage type	Accessories/DL
EAN code (Model)	4045761190450
Predecessor Model	ODC 824
Customs tariff number	85258900
CE Label	yes

Technical data - Primary

Camera type	USB2 Microscope camera - C-Mount
Camera resolution	5,1 MP
Sensor type	Aptina CMOS

Sensor size 1/2,5"

Technical data - Construction

Color capability	Colour
Camera mounting	C-Mount
Interfaces	USB 2.0
Operating system	Win 10;Win 11;Win 7;Win 8
Applicable software	Microscope VIS 2.0 Pro (OXM 902)
Software functions	Bildaufnahme;Bildbearbeitung/Dokumentation;Messungen;Videoaufnahme
Dimensions (W×D×H)	80×80×45 mm

Technical data - Power supply

Power Plug	USB B
Protection Class	Class III - Low voltage

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	164×164×75 mm
Gross weight	1.00 kg
Net weight	0.60 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	342.022 g
Packaging component - by weight - plastics	12.0959 g
Packaging component - by weight - styrofoam	0 g
Packaging component - by weight - foam	62.9821 g
Delivery time	1 d

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Verification

WEEE	yes
------	-----

FUNCTIONS

Standard



ACCESSORIES

Model	Description
ODC-A8102	Eyepiece Adapter KERN ODC-A8102

Model	Description
ODC-A8103	Eyepiece Adapter KERN ODC-A8103
ODC-A8105	Eyepiece Adapter KERN ODC-A8105
OXM 902	Camera software VIS KERN microscope OXM 902

SPARE PARTS

Model	Description
ODC-A8201	USB 2.0 cable ODC-A8201

RELATED ARTICLES DATA

Model name	Camera type	Camera resolution	Sensor type	Sensor size
ODC 825	USB2 Microscope camera - C-Mount	5,1 MP	Aptina CMOS	1/2,5"