

PRODUCT BROCHURE

KERN SET-311

Date: 18.03.2025

ENGLISH

Software EasyTouch SET-311 **Quick Key**






EasyTouch Quick Key – select master data objects quickly and user-related

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

PRODUCT SHOWCASE



DESCRIPTION

- Note: Prerequisite for this set is the basic program SET-01 Base
- Quick Keys can be created in numerous sets, e. g.
 - SET-01 Base
 - SET-05 Dynamic
 - SET-06 Classification
 - SET-07 Variable
 - SET-08 Take Out
 - SET-31 Tolerance
 - SET-32 Count
 - SET-33 Target Count
- Quick selection keys for master data: Through ET Quick Keys, specific master data objects can be permanently displayed for quick selection. The Quick Keys are displayed on the screen immediately
- Saving: Whether selected in the conventional way or by using Quick Keys, every result of a weighing or piece count will be stored in relation to the selected master data object
- Selecting master data objects: The specific master data objects, which are displayed as Quick Keys, can either

be selected individually or as a defined master data group, which contains the master data as Quick Keys

- User-related: The Quick Keys can be created as a scheme which relates to particular user roles. This means that the users of a particular role are shown the relevant master data objects for quick selection. This allows, for example, different departments in the same company to work with different master data without overlapping
- Saving and printing: For each scheme you can set whether, in addition to the being saved, it should also be automatically printed out or not
- Communication between balance/terminal device: The balance(s) can communicate with the PC, laptop or tablet by serial connection, USB, Bluetooth, Ethernet or WiFi

License information

- User: You can store as many users as you need in one license
- Balances: You can store and operate as many balances as you need in one license
- Communication between balance/terminal device: The balance(s) can communicate with the PC, laptop or tablet by serial connection, USB, Bluetooth, Ethernet or WiFi

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	SET-311
Model-Serie	ETS
Assortment	KERN
Product Type	Software
Product group	EasyTouch
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1486546
Customs tariff number	85234910

Standard



Option



Technical data - Packing & Shipping

Delivery time	1 d
---------------	-----

FUNCTIONS