

## PRODUCT BROCHURE

---

# SAUTER AD 003 1

Date: 30.04.2025

ENGLISH



Image not  
available

Screw-in tension clamp for AD 0031, Aluminium, = 2 kN, 1 piece included in the scope of delivery

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                                                        |
|-----------------------|--------------------------------------------------------|
| Model article-code    | AD 0031                                                |
| Model-Serie           | AD                                                     |
| Assortment            | Sauter                                                 |
| Product Type          | Force measurement                                      |
| Product group         | Clamps & clamping tools                                |
| Product usage type    | Trading parts list<br>element;Accessories/Se<br>rvices |
| EAN code (Model)      | 4045761174450                                          |
| Customs tariff number | 90278990                                               |

|                                 |              |
|---------------------------------|--------------|
| Dimensions packaging<br>(W×D×H) | 120×75×30 mm |
|---------------------------------|--------------|

|              |         |
|--------------|---------|
| Gross weight | 0.70 kg |
|--------------|---------|

|            |         |
|------------|---------|
| Net weight | 0.70 kg |
|------------|---------|

|                 |                |
|-----------------|----------------|
| Shipping method | Parcel service |
|-----------------|----------------|

|           |                    |
|-----------|--------------------|
| UN number | No dangerous goods |
|-----------|--------------------|

|                                               |     |
|-----------------------------------------------|-----|
| Packaging component -<br>by weight - plastics | 0 g |
|-----------------------------------------------|-----|

|               |      |
|---------------|------|
| Delivery time | 14 d |
|---------------|------|

ENGLISH

### Technical data - Primary

|                  |        |
|------------------|--------|
| Max. Force       | 2 kN   |
| Max. Clamp width | 100 mm |

### Technical data - Construction

|                              |              |
|------------------------------|--------------|
| Material                     | aluminium    |
| Dimension housing<br>(W×D×H) | 120×75×30 mm |

### Technical data - Environmental conditions

|                                |        |
|--------------------------------|--------|
| Maximum storage<br>temperature | 40 °C  |
| Minimum storage<br>temperature | -10 °C |

### Technical data - Packing & Shipping

## ACCESSORIES

| Model   | Description       |
|---------|-------------------|
| AD 0034 | Jaws KERN AD 0034 |
| AD 0035 | Jaws KERN AD 0035 |
| AD 0036 | Jaws KERN AD 0036 |
| AD 0037 | Jaws KERN AD 0037 |
| AD 0038 | Jaws KERN AD 0038 |
| AD 0039 | Jaws KERN AD 0039 |
| AD 0040 | Jaws KERN AD 0040 |
| AD 0041 | Jaws KERN AD 0041 |
| AD 0042 | Jaws KERN AD 0042 |
| AD 0043 | Jaws KERN AD 0043 |
| AD 0044 | Jaws KERN AD 0044 |
| AD 0045 | Jaws KERN AD 0045 |
| AD 0046 | Jaws KERN AD 0046 |