

KERN ORA 6HA

Refractive index measurement for laboratories and the industry



Article codes, product groups & Status information

Model article-code	ORA 6HA
Model-Serie	ORA-H
Assortment	Optics
Product Type	Refractometers
Product group	Analogue refractometers
Product usage type	Main article
EAN code (Model)	4045761483491
Predecessor Model	ORA 6HA-2023a; ORA 6HB
Planned start date	2025-1HJ
Customs tariff number	90275000

Technical data - Primary

Scales	Water content
Scale Range	12 % - 30 %
ATC	Yes

Technical data - Construction

Type	Handheld refractometer
Display scale	Analogue
Dimensions (W×D×H)	168×40×40 mm

Technical data - Measuring system

Measuring method	Light refraction
Measuring Temp min	10 °C
Measuring Temp max	30 °C
Readout	0.1 %

Functions

Standard

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	215×94×65 mm
Gross weight	0.45 kg
Net weight	0.30 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	26.14 g
Packaging component - by weight - plastics	26.14 g
Delivery time	1 d

KERN ORA 6HA

Refractive index measurement for laboratories and the industry

