

PRODUCT BROCHURE

KERN 323-024

Date: 09.03.2026

ENGLISH



Weight set, F1, 1 mg - 50 g, Knob, stainless steel polished, in plastic case

CONTACT


KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



DESCRIPTION

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	323-024
Model-Serie	323
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main article
EAN code (Model)	404576 108075 1
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 mg - 50 g
---------------	-------------

Case plastic	313-052-400
Weight set case	in plastic case

Technical data - Construction

Material	stainless steel polished
Is weight set	yes
Number of weights in set	20
Type of case	in plastic case

Technical data - Measuring system

Nominal value	1 mg - 50 g
OIML class	F1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	250×210×75 mm
Gross weight	0.50 kg
Net weight	0.50 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-401
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>962-401</u>	DAkkS Calibration KERN 962-401
<u>962-401R</u>	DAkkS Recalibration KERN 962-401R
<u>921-401</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-401

Model	Description
<u>920-702</u>	Marking (Last digits of serial number) KERN 920-702

SPARE PARTS

Model	Description
<u>313-052-400</u>	Plastic Weight Box KERN 313-052-400