


PRODUCT BROCHURE

SAUTER AC 18R

Date: 30.04.2025

ENGLISH



Screw tension clamp with 2 wide jaw grips for tension and extraction tests

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	AC 18R
Model-Serie	AC
Assortment	Sauter
Product Type	Force measurement
Product group	Clamps & clamping tools
Product usage type	Trading parts list element;Accessories/Services
EAN code (Model)	4045761122796
Customs tariff number	90278990

Technical data - Primary

Max. Force	5 kN
Max. Clamp width	60 mm

Technical data - Construction

Material	steel
----------	-------

Dimension housing (W×D×H) 160×100×140 mm

Technical data - Measuring system

Mounting threads External thread M10

Technical data - Packing & Shipping

Dimensions packaging (W×D×H) 160×100×140 mm

Gross weight 0.80 kg

Net weight 0.75 kg

Shipping method Parcel service

UN number No dangerous goods

Packaging component - by weight - plastics 10 g

Delivery time 1 d

FUNCTIONS

Standard



ACCESSORIES

Model	Description
AFM 03	Sleeve SAUTER AFM 03
AFM 06	Adapter KERN AFM 06
AFM 11	assembly and repair AFM 11
AFM 16	Adapter KERN AFM 16