

KERN 366-94

Individual weight, OIML Class M3, 1 kg, Knob, cast iron lacquered



Article codes, product groups & Status information

Model article-code	366-94
Model-Serie	366
Assortment	KERN
Product Type	Test weights
Product group	OIML M3 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761482081
Predecessor Model	366-94-2024a
Planned start date	2025-2HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	1 kg
Tolerance (OIML)	500 mg

Technical data - Construction

Material	Cast iron lacquered
----------	---------------------

Technical data - Measuring system

Nominal value	1 kg
OIML class	M3
OIML conform	Yes
Tolerance (OIML)	500 mg
Density (kg/m ³)	7100
Volume (cm ³)	140,8451

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	55×55×87 mm
Gross weight	1.00 kg
Net weight	1.00 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	No
-----------------------	----

Optional Services

DAkKS Calibration	962-641
-------------------	---------

Functions

KERN 366-94

Individual weight, OIML Class M3, 1 kg, Knob, cast iron lacquered



Standard

