

KERN OZP 551

Individuality, variety and flexible working through our modular construction system ?
Stereo microscope heads



Article codes, product groups & Status information

Model article-code	OZP 551
Model-Serie	OZB-M
Assortment	Microscopes & Optical Metrology
Product Type	Microscopes
Product group	Stereo microscopes
Product usage type	Spare part;Main article
EAN code (Model)	4045761183643
Successor Model	OZP 151
Planned start date	2015-2HJ
Customs tariff number	90119000
Model Status	Yes
Country of Origin	CN
Selective distribution	0

Technical data - Primary

Optical System	Greenough
Tube	Binocular
Eyepiece	Eyepiece HSWF 10x / Ø 23 mm with anti-fungus, high eye point
Objectives	0,6x - 5,5x

Technical data - Construction

Tube angle	35°
Tube rotation	Yes
Eye distance (max)	76 mm
Packaging component - by weight - foam	0 g
Delivery time	1 d

Eye distance (min)	52 mm
Adjustable dioptre	Both-sided
Visual field min	38,3 mm
Visual field max	4,2 mm
Working distance	108 mm
Sample height max	110 mm
Max resolution	0,0042 mm

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	310×300×220 mm
Gross weight	3,0 kg
Net weight	2,6 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	262,779 g
Packaging component - by weight - plastics	8,7776 g
Packaging component - by weight - styrofoam	277,043 g

KERN OZP 551

Individuality, variety and flexible working through our modular construction system ?
Stereo microscope heads



Functions

Standard

