

KERN 317-08

Individual weight, OIML Class E2, 200 g, Knob, stainless steel polished



Article codes, product groups & Status information

Model article-code	317-08
Model-Serie	317
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761312562
Predecessor Model	317-8-2.021a; 317-8-2.022i
Planned start date	2022-2HJ
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	200 g
Tolerance (OIML)	0,3 mg
Case aluminium	317-080-600
Case plastic	317-080-400
Case wood	317-080-100

Technical data - Construction

Material	Stainless steel polished
----------	--------------------------

Technical data - Measuring system

Nominal value	200 g
OIML class	E2
OIML conform	Yes

Functions

Standard

Tolerance (OIML) 0,3 mg

Density (kg/m³) 7.950

Volume (cm³) 25

Technical data - Packing & Shipping

Dimensions packaging (W×D×H) 28×28×47 mm

Gross weight 0,20 kg

Net weight 0,20 kg

Shipping method Parcel service

Packaging component - by weight - plastics 0 g

Delivery time 1 d

Optional Services

DAkKS Calibration 962-338

KERN 317-08

Individual weight, OIML Class E2, 200 g, Knob, stainless steel polished



Option

