KERN MWS 300K-1LM

Stretcher scale with low overall height for easy access



Article codes, product groups & Status information	
Model article-code	MWS 300K-1LM
Model-Serie	MWS
Assortment	Scales
Segment	Medical products
Product group	Wheelchair platform scales
Product usage type	Main article
EAN code (Model)	4045761192881
Predecessor Model	MWS 300K-1LM-2.016a
Successor Model	MWN 300K-1LM
Planned start date	2017-2HJ
Customs tariff number	84231090
Model Status	Yes
CE Label	Yes
Country of Origin	CN

Technical data - Primary	
Weighing capacity	300 kg
Readability	0,1 kg
Verification value (e)	0,1 kg

Technical data - Construction	
Material weighing plate	Steel, lacquered
Bubble Level	Yes
Revolving screw feet	Yes
Weighing surface (d)	800×1.200 mm
Dimensions of weighing plate (W×D)	800×1.200 mm
Dimension completely mounted (W×D×H)	1.506×920×68 mm
Material casing	Steel, lacquered

Technical data - Display	
Display height	25 mm

Technical data - Power supply	
Battery	6×15 V AA
Operating time	50 h
Mains adapter/adapter included	CH;EURO;UK
Charging time	14 h
Input Voltage	100 V - 240 V, 50 / 60 Hz
Operating time (Backlight off)	45 h
Rechargeable battery	Rchrg. battery optional - intern
Interfaces	RS-232 standard

Technical data - Measuring system	
Calibration method	Adjusting with external weight
Linearity	± 0,1 kg
Readability	0,1 kg
Recommended adjusting weight	300 kg (M1)
Repeatability	0,1 kg
Weighing units	Kg
Weighing system	Strain gauge
Max. divergency at eccentric load (1/3 max)	0,1 kg

KERN MWS 300K-1LM

Stretcher scale with low overall height for easy access

Technical data - Environmental conditions	
Maximum humidity	80 %
Maximum operating temperature	35 °C
Minimum ambient temperature	5 °C

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	1.650×1.008×123 mm
Gross weight	55 kg
Net weight	42 kg
Shipping method	Freight
Packaging component - by weight - cardboard	6,5 kg
Packaging component - by weight - styrofoam	6,5 kg
Delivery time	2 d

Technical data - Verification	
Verification Value E	0,1 kg
Minimum load (Min)	2 kg
Verification class	III
MPG class	1m
CND code	V0399
GMDN code	15053

Optional Services	
DAkkS Calibration	963-129

Functions

Standard



Option

