KERN 326-08

Individual weight, OIML Class F1, 200 g, Eco Shape, stainless steel polished



Article codes, product groups & Status information	
Model article-code	326-08
Model-Serie	326
Assortment	Scales
Segment	Test weights
Product group	OIML F1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761247291
Predecessor Model	326-8-2.018a; 326-8-2.025a; 327-68; 329-8
Planned start date	2026-2HJ
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary	
Nominal value	200 g
Tolerance (OIML)	1 mg
Case aluminium	317-080-600
Case plastic	347-080-400
Case wood	317-080-100

Technical data - Construction	
Material	Stainless steel polished

Technical data - Measuring system	
Nominal value	200 g
OIML class	F1
OIML conform	Yes
Tolerance (OIML)	1 mg
Density (kg/m³)	7.960
Volume (cm³)	25

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	28×28×47 mm
Gross weight	0,20 kg
Net weight	0,20 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services	
DAkkS Calibration	962-438

Functions

Standard

KERN 326-08

Individual weight, OIML Class F1, 200 g, Eco Shape, stainless steel polished



Option

