KERN 347-095-N00

Individual weight (Newton), 5 N, Slot weight (Type A), stainless steel fine turned



Article codes, product groups & Status information	
Model article-code	347-095-N00
Model-Serie	347-5-N
Assortment	Scales
Segment	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761277281
Planned start date	2022-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary	
Tolerance (OIML)	25 mg
Case plastic	347-110-400

Technical data - Construction	
Material	Stainless steel fine turned

Technical data - Measuring system	
OIML class	M1
OIML conform	No
Density (kg/m^3)	7.950
Volume (cm^3)	59,5238

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	48×48×40 mm
Gross weight	0,50 kg
Net weight	0,50 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	6 d

Optional Services	
DAkkS Calibration	962-639N

Functions

KERN 347-095-N00

Individual weight (Newton), 5 N, Slot weight (Type A), stainless steel fine turned

Standard



Option

