KERN MPK 200K-1P

Personal scale with column for convenient, ergonomically-optimised operation



Article codes, product groups & Status information	
Model article-code	MPK 200K-1P
Model-Serie	MPK
Assortment	Scales
Segment	Medical products
Product group	Bathroom scales
Product usage type	Main article
EAN code (Model)	4045761201989
Planned start date	2018-1HJ
Customs tariff number	84231090
CE Label	Yes
Country of Origin	CN

Technical data - Primary	
Weighing capacity	250 kg
Readability	0,1 kg
Dimension weighing platform (W×D×H)	350×350×45 mm

Technical data - Construction	
Material weighing plate	Plastic
Weighing surface (d)	350×350 mm
Dimensions of weighing plate (W×D)	350×350 mm
Dimension completely mounted (W×D×H)	350×420×955 mm

Technical data - Display	
Display height	28 mm

Technical data - Power supply	
Battery	4×15 V AA
Operating time	100 h

Technical data - Measuring system	
Calibration method	Adjusting with external weight
Linearity	± 1 kg
Readability	0,1 kg
Recommended adjusting weight	250 kg (M1)
Repeatability	0,1 kg
Weighing units	Kg;lb
Weighing system	Strain gauge

Technical data - Environmental conditions	
Maximum humidity	80 %
Maximum operating temperature	40 °C
Minimum ambient temperature	10 °C
Minimum humidity	5 %

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	445×1.003×145 mm
Gross weight	7 kg
Net weight	3,8 kg
Shipping method	Parcel service

Packaging component - by weight - paper	604 g
Packaging component - by weight - cardboard	1.097 g
Packaging component - by weight - plastics	17 g
Packaging component - by weight - styrofoam	394 g
Packaging component - by weight - wood	370 g
Delivery time	1 d

Optional Services

DAkkS Calibration	963-129

Functions

Standard

