KERN MGC 100K-1

The robust classic for quick determination of body weight



Article codes, product groups & Status information	
Model article-code	MGC 100K-1
Model-Serie	MGC
Assortment	KERN
Product Type	Medical products
Product group	Bathroom scales
Product usage type	Spare part;Trading parts list element;Main item
EAN code (Model)	4045761202078
Planned start date	2018-1HJ
Customs tariff number	84231090
Country of Origin	CN

Technical data - Primary	
Weighing capacity	150 kg
Readability	1 kg
Dimension weighing platform (W×D×H)	265×210×40 mm

Technical data - Construction	
Material weighing plate	Plastic
Weighing surface (d)	265×210 mm
Dimensions of weighing plate (W×D)	265×210 mm
Dimension completely mounted (W×D×H)	290×425×90 mm

Technical data - Display	
Display height	8 mm

Technical data - Measuring system	
Calibration method	Calibration not possible
Linearity	± 2 kg
Readability	1 kg
Repeatability	1 kg
Weighing units	Kg
Weighing system	Mechanische Waage

Technical data - Environmental conditions	
Maximum humidity	80 %
Maximum operating temperature	40 °C
Minimum ambient temperature	-10 °C
Minimum humidity	5 %

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	470×315×113 mm
Gross weight	3,6 kg
Net weight	3,2 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	300 g
Delivery time	1 d

Functions

Standard

KERN MGC 100K-1

The robust classic for quick determination of body weight

