KERN MBC 15K2DM

Elegantly-shaped baby scale – with approval for professional mobile or stationary medical use in medical diagnostics, verification optional



Article codes, product groups & Status information	
Model article-code	MBC 15K2DM
Model-Serie	MBC
Assortment	Scales
Segment	Medical products
Product group	Baby balances
Product usage type	Main article
EAN code (Model)	4045761192935
Predecessor Model	MBC 15K2DM-2.017a
Planned start date	2017-2HJ
Customs tariff number	84231090
CE Label	Yes
Country of Origin	CN

Technical data - Primary	
Weighing capacity	6 kg;15 kg
Readability	0,002 kg;0,005 kg
Verification value (e)	0,002 kg;0,005 kg

Technical data - Construction	
Material weighing plate	Plastic
Bubble Level	Yes
Weighing surface (d)	600×276 mm
Dimensions of weighing plate (W×D)	600×276 mm
Dimension completely mounted (W×D×H)	600×410×112 mm
Material casing	Plastic

Technical data - Display	
Display height	25 mm

Technical data - Power supply	
Battery	6×15 V AA
Operating time	20 h
Mains adapter/adapter included	CH;EURO;UK
Charging time	12 h
Input Voltage	100 V - 240 V, 50 / 60 Hz
Operating time (Backlight off)	40 h
Rechargeable battery	Rchrg. battery optional - intern

Technical data - Measuring system	
Calibration method	Adjusting with external weight
Linearity	± 0,002 kg;0,005 kg
Recommended adjusting weight	15 kg (M1)
Repeatability	0,002 kg;0,005 kg
Weighing units	Kg
Weighing system	Strain gauge

Technical data - Environmental conditions	
Maximum humidity	80 %
Maximum operating temperature	40 °C
Minimum ambient temperature	-10 °C

KERN MBC 15K2DM

Elegantly-shaped baby scale – with approval for professional mobile or stationary medical use in medical diagnostics, verification optional


Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	480×654×177 mm
Gross weight	6 kg
Net weight	4,6 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	301 g
Packaging component - by weight - styrofoam	1,204 g
Delivery time	1 d

Technical data - Verification	
Minimum load (Min)	0,04 kg;0,1 kg
Verification by KERN possible	Yes
Verification class	III
MPG class	1m
CND code	Z12080402
GMDN code	35324

Optional Services	
DAkkS Calibration	963-128
Verification	965-128.

Functions

Standard



Option



Factory

