PRODUCT BROCHURE

KERN RH-A01

Date: 30.05.2026

ENGLISH



Disposable sample plates aluminium, Ø 92 mm, unit of 50 pieces

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

Scan here for more information

DESCRIPTION

- Suitable for moisture analyzer KERN DAB, DBS, DIS, DLB

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	RH-A01
Model-Serie	RH
Assortment	KERN
Product Type	Laboratory balances
Product group	Accessories (Printer, etc)
Product usage type	Accessories/DL
EAN code (Model)	4045761503588
Predecessor Model	MLB-A01A
Customs tariff number	90160090
Country of Origin	CN

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	105×105×80 mm
Gross weight	0.35 kg
Net weight	0.30 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	37 g
Packaging component - by weight - plastics	10 g
Delivery time	1 d

FUNCTIONS

Standard



ENGLISH