
KERN MWC 600K-1

Mobile wheelchair platform scale with handles and rollers



Article codes, product groups & Status information	
Model article-code	MWC 600K-1
Model-Serie	MWC
Assortment	Scales
Segment	Medical products
Product group	Wheelchair platform scales
Product usage type	Main article
EAN code (Model)	4045761482869
Predecessor Model	MWB 300K-1
Planned start date	2025-1HJ
Customs tariff number	84231090
CE Label	Yes
Country of Origin	CN

Technical data - Primary	
Weighing capacity	600 kg
Readability	0,1 kg

Technical data - Construction	
Material weighing plate	Steel, powder-coated
Bubble Level	Yes
Revolving screw feet	Yes
Weighing surface (d)	800×800 mm
Dimensions of weighing plate (W×D)	800×800 mm
Dimension completely mounted (W×D×H)	1.100×1.000×80 mm
Material casing	Steel, powder-coated
IP protection indicator	No IP protection
IP protection platform	IP67
IP protection load cell	IP67

Technical data - Display	
Display height	25 mm

Technical data - Power supply	
Battery	4×15 V AA
Operating time	60 h
Mains adapter/adapter included	AUS;CH;EURO;UK;US
Input Voltage	100 V - 240V AC 50/60Hz 0,3A

Technical data - Measuring system	
Calibration method	Adjusting with external weight
Linearity	± 0,6 kg
Readability	0,1 kg
Recommended adjusting weight	600 kg (M1)
Repeatability	0,2 kg
Weighing units	Kg;lb
Weighing system	Strain gauge
Max. divergency at eccentric load (1/3 max)	1 kg

Technical data - Environmental conditions	
Maximum humidity	80 %

Maximum operating temperature	40 °C
Minimum ambient temperature	-10 °C
Minimum humidity	5 %

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	1.135×1.050×200 mm
Gross weight	65 kg
Net weight	40 kg
Shipping method	Freight
Packaging component - by weight - wood	22 kg
Delivery time	2 d

Functions

Standard

