





Software EasyTouch
SET-21
Formulation

Upgrade your laboratory balance to a high-end recipe balance for the pharmaceutical, chemical and food industries

Features

- By combining the KERN PCB with the Easy-Touch SET-21 you can manage recipes and save them against a user role and in a way that they cannot be altered. When implementing a recipe, EasyTouch Formulation guides the user through the recipe, component by component and provides support for accurate dispensing through the optical and acoustic dosing assistant
- Preventing errors: EasyTouch helps you to avoid errors such as incorrect components or inaccurate weight information to a large extent
- The correction function recognises overfilling and indicates this. If it is no longer possible to perform a take-out, then Easytouch can automatically adjust the recipe and this prevents losses of expensive or limited resources
- The bundle contains:
  - Precision balance KERN
    1 PCB 300-3 / PCB 1000-2
    2 PCB 10000-1
  - Software EasyTouch SET-01 Base
  - Software EasyTouch SET-21 Formulation
  - USB interface adapter KERN KUP-03

STANDARD

OPTION

Model	Weighing capacity [Max]	Readability [d]	Reproducibility	Resolution	Weighing surface	Overall dimensions W×D×H mm	Net weight approx.	Options Calibr. Certificate DAkkS accr. KERN
KERN	g	g	g	Points	mm		kg	
PCB 300-3-REZ	360	0,001	0,002	360.000	105	163×245×65	0,85	963-127
PCB 1000-2-REZ	1200	0,01	0,01	120.000	130×130	163×245×65	1,2	963-127
PCB 10000-1-REZ	10000	0,1	0,1	100.000	150×170	163×245×65	1,4	963-128