KERN OZL 461

Individuality, variety and flexible working through our modular construction system ?
Stereo microscope heads



Article codes, product groups & Status information	
Model article-code	OZL 461
Model-Serie	OZB-M
Assortment	Optics
Product Type	Microscopes
Product group	Stereo microscopes
Product usage type	Main article
EAN code (Model)	4045761189348
Planned start date	2016-1HJ
Customs tariff number	90119000
CE Label	Yes
Country of Origin	CN

Technical data - Primary	
Optical System	Greenough
Tube	Binocular
Eyepiece	Eyepiece HWF 10x / Ø 20mm High Eye Point
Objectives	0,7x - 4,5x

Technical data - Construction	
Tube angle	45°
Tube rotation	Yes
Eye distance (max)	75 mm
Eye distance (min)	55 mm
Adjustable dioptre	Both-sided
Packaging component - by weight - foam	0 g
Delivery time	1 d

Technical data - Verification	
WEEE	Yes
Visual field min	28,6 mm
Visual field max	4,4 mm
Working distance	105 mm
Sample height max	140 mm
Max resolution	3,7 mm

Technical data - Power supply	
Protection Class	No protection class

Technical data - Environmental conditions	
Storage Temp Min	-5 °C
Stotage Temp Max	40 °C

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	310×245×230 mm
Gross weight	2,0 kg
Net weight	1,2 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	799,71 g
Packaging component - by weight - plastics	14,64 g
Packaging component - by weight - styrofoam	99,735 g

Functions

Standard

