KERN SFE 300K-1LNM

Platform scale with IP65 dust and spray protection, verification optional



Article codes, product groups & Status information	
Model article-code	SFE 300K-1LNM
Model-Serie	SFE
Assortment	Scales
Segment	Industrial balance
Product group	Platform scales
Product usage type	Main article
EAN code (Model)	4045761244726
Predecessor Model	SFE 300K-1LM; SFE 300K-1LNM-2.021a
Planned start date	2021-1HJ
Customs tariff number	84238220
CE Label	Yes
Country of Origin	CN

Technical data - Primary	
Weighing capacity	300 kg
Readability	0,1 kg
Net weight	20,2 kg
Platform size	650×500 mm
Indicator cable length	3 m
Dimension weighing platform (W×D×H)	650×500×139 mm

Technical data - Construction	
Dimensions display casing (W×D×H)	195×118×83 mm
Material weighing plate	Stainless steel
Material Platform construction	Steel, powder-coated
Bubble Level	Yes
Revolving screw feet	Yes
Weighing surface (d)	650×500 mm
Dimensions of weighing plate (W×D)	650×500 mm
Indicator cable length	3 m
Wall mount	Yes
IP protection indicator	IP65
IP protection platform	No IP protection
IP protection load cell	IP65

Technical data - Display	
Display height	22 mm
Backlight	Yes

Technical data - Power supply	
Mains adapter/adapter included	CH;EURO;UK;US
Charging time	12 h
Input Voltage	220 V - 240 V AC, 50 Hz
Operating time (Backlight off)	35 h

Technical data - Measuring system	
Calibration method	Adjusting with external weight
Linearity	± 0,1 kg
Readability	0,1 kg

KERN SFE 300K-1LNM

Platform scale with IP65 dust and spray protection, verification optional

Recommended adjusting weight	300 kg (M1)
Repeatability	0,1 kg
Resolution	3.000
Stabilisation time under laboratory conditions	2 s
Tare range	299,9 kg
Warm up time	10 min
Weighing system	Strain gauge
Weighing units	Kg
load cell connection	6-Leiter

Technical data - Environmental conditions

Maximum humidity	80 %
Maximum operating temperature	40 °C
Minimum ambient temperature	-10 °C
Minimum humidity	0 %

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	1.050×675×250 mm
Gross weight	24 kg
Net weight	22 kg
Shipping method	Parcel service
UN number	3496
Packaging component - by weight - cardboard	1,05 g
Packaging component - by weight - styrofoam	1,05 g
Delivery time	1 d

Technical data - Verification

Verification value (e)	0,1 kg
Minimum load (Min)	2 kg
Verification by KERN possible	Yes
Verification class	III

Optional Services

DAkkS Calibration	963-129
Verification	965-229

Functions

Standard



Option



Factory

