KERN HDB 10K10N

Practical hanging scale for rapid, mobile weighing



Article codes, product groups & Status information	
Model article-code	HDB 10K10N
Model-Serie	HDB
Assortment	Scales
Product Type	Industrial balance
Product group	Hanging balances
Product usage type	Main article
EAN code (Model)	4045761135536
Planned start date	2011-1HJ
Customs tariff number	84238129
CE Label	Yes
Country of Origin	CN

Technical data - Primary	
Weighing capacity	10 kg
Readability	0,01 kg
Dimensions casing (W×D×H)	72×25×105 mm
Net weight	0,189 kg

Technical data - Construction	
Material casing	Plastic
Dimensions casing	72×25×105 mm
Internal measurement frequency	5 Hz
Material load-hook	Steel
TUEV (TÜV)	Yes
IP protection load cell	IP65

Technical data - Display	
Display height	12 mm
Backlight	No

Technical data - Power supply	
Battery	2×15 V AAA
Operating time	100 h
Operating time (Backlight off)	100 h

Technical data - Measuring system	
Calibration method	Adjusting with external weight
Linearity	± 0,03 kg
Readability	0,01 kg
Recommended adjusting weight	10 kg (M1)
Repeatability	0,01 kg
Resolution	1.000
Warm up time	10 min
Weighing units	Kg;lb;N

Technical data - Environmental conditions	
Maximum humidity	90 %
Maximum operating temperature	35 °C
Minimum ambient temperature	5 °C

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	225×150×70 mm

KERN HDB 10K10N


Practical hanging scale for rapid, mobile weighing

Gross weight	0,45 kg
Net weight	0,20 kg
Shipping method	Parcel service
Packaging component - by weight - paper	17 g
Packaging component - by weight - cardboard	690 g
Packaging component - by weight - plastics	1 g
Packaging component - by weight - styrofoam	77 g
Packaging component - by weight - wood	16 g
Delivery time	1 d

Optional Services	
DAkkS Calibration	963-128H

Functions

Standard



Option

